

RPB7 Polyclonal Antibody

Cat No: HR1AP10717

For research use only

Overview

Product Name	RPB7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	POLR2G RPB7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RNA polymerase II subunit G(POLR2G) Homo sapiens This gene encodes the seventh largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The protein functions in transcription initiation, and is also thought to help stabilize transcribing polymerase molecules during elongation. [provided by RefSeq, Jan 2009],