

RPB3 Polyclonal Antibody

Cat No: HR1AP10715

For research use only

Overview

Product Name	RPB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	POLR2C A-152E5.7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RNA polymerase II subunit C(POLR2C) Homo sapiens This gene encodes the third largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene contains a cysteine rich region and exists as a heterodimer with another polymerase subunit, POLR2J. These two subunits form a core subassembly unit of the polymerase. A pseudogene has been identified on chromosome 21. [provided by RefSeq, Jul 2008],</p>