

RPA1 Polyclonal Antibody

Cat No: HR1AP10712

For research use only

Overview

Product Name	RPA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	189kDa
GeneID?Human?	POLR1A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RNA polymerase I subunit A(POLR1A) Homo sapiens The protein encoded by this gene is the largest subunit of the RNA polymerase I complex. The encoded protein represents the catalytic subunit of the complex, which transcribes DNA into ribosomal RNA precursors. Defects in this gene are a cause of the Cincinnati type of acrofacial dysostosis. [provided by RefSeq, May 2016],