

SMBP2 Polyclonal Antibody

Cat No: HR1AP10711

For research use only

Overview

Product Name	SMBP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	109kDa
GeneID?Human?	IGHMBP2 SMBP2 SMUBP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	immunoglobulin mu binding protein 2(IGHMBP2) Homo sapiens This gene encodes a helicase superfamily member that binds a specific DNA sequence from the immunoglobulin mu chain switch region. Mutations in this gene lead to spinal muscle atrophy with respiratory distress type 1. [provided by RefSeq, Jul 2008],