

DNA2 Polyclonal Antibody

Cat No: HR1AP10709

For research use only

Overview

Product Name	DNA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	116kDa
GeneID?Human?	DNA2 DNA2L KIAA0083
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>DNA replication helicase/nuclease 2(DNA2) Homo sapiens This gene encodes a member of the DNA2/NAM7 helicase family. The encoded protein is a conserved helicase/nuclease involved in the maintenance of mitochondrial and nuclear DNA stability. Mutations in this gene are associated with autosomal dominant progressive external ophthalmoplegia-6 (PEOA6) and Seckel syndrome 8. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2014],</p>