

## CDT1 Polyclonal Antibody

Cat No: HR1AP10707

For research use only

### Overview

Product Name	CDT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	CDT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromatin licensing and DNA replication factor 1(CDT1) Homo sapiens The protein encoded by this gene is involved in the formation of the pre-replication complex that is necessary for DNA replication. The encoded protein can bind geminin, which prevents replication and may function to prevent this protein from initiating replication at inappropriate origins. Phosphorylation of this protein by cyclin A-dependent kinases results in degradation of the protein. [provided by RefSeq, Mar 2011],</p>