

## PSF2 Polyclonal Antibody

Cat No: HR1AP10706

For research use only

### Overview

Product Name	PSF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	GINS2 PSF2 CGI-122 DC5 HSPC037
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GINS complex subunit 2(GINS2) Homo sapiens The yeast heterotetrameric GINS complex is made up of Sld5 (GINS4; MIM 610611), Psf1 (GINS1; MIM 610608), Psf2, and Psf3 (GINS3; MIM 610610). The formation of this complex is essential for the initiation of DNA replication in yeast and Xenopus egg extracts (Ueno et al., 2005 [PubMed 16287864]). See GINS1 for additional information about the GINS complex.[supplied by OMIM, Mar 2008],