

DPOA2 Polyclonal Antibody

Cat No: HR1AP10702

For research use only

Overview

Product Name	DPOA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	POLA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:The N-terminal 240 amino acids are sufficient to mediate complex formation.,function:May play an essential role at the early stage of chromosomal DNA replication by coupling the polymerase alpha/primase complex to the cellular replication machinery.,PTM:Phosphorylated in a cell cycle-dependent manner, in G2/M phase.,similarity:Belongs to the DNA polymerase alpha subunit B family.,subunit:DNA polymerase alpha:primase is a four subunit enzyme complex, which is assembled throughout the cell cycle, and consists of the two DNA polymerase subunits A and B, and the DNA primase large and small subunits. Subunit B binds to subunit A.,</p>