

ABC3H Polyclonal Antibody

Cat No: HR1AP10701

For research use only

Overview

Product Name	ABC3H Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	APOBEC3H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>apolipoprotein B mRNA editing enzyme catalytic subunit 3H(APOBEC3H) Homo sapiens This gene encodes a member of the apolipoprotein B mRNA-editing enzyme catalytic polypeptide 3 family of proteins. The encoded protein is a cytidine deaminase that has antiretroviral activity by generating lethal hypermutations in viral genomes. Polymorphisms and alternative splicing in this gene influence its antiretroviral activity and are associated with increased resistance to human immunodeficiency virus type 1 infection in certain populations. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009],</p>