

## DLG3 Polyclonal Antibody

Cat No: HR1AP10700

For research use only

### Overview

Product Name	DLG3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	DLG3 KIAA1232
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	discs large MAGUK scaffold protein 3(DLG3) Homo sapiens This gene encodes a member of the membrane-associated guanylate kinase protein family. The encoded protein may play a role in clustering of NMDA receptors at excitatory synapses. It may also negatively regulate cell proliferation through interaction with the C-terminal region of the adenomatosis polyposis coli tumor suppressor protein. Mutations in this gene have been associated with X-linked mental retardation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Oct 2009],