

DGAT2 Polyclonal Antibody

Cat No: HR1AP10699

For research use only

Overview

Product Name	DGAT2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	DGAT2 HMFN1045 UNQ738/PRO1433
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>diacylglycerol O-acyltransferase 2(DGAT2) Homo sapiens This gene encodes one of two enzymes which catalyzes the final reaction in the synthesis of triglycerides in which diacylglycerol is covalently bound to long chain fatty acyl-CoAs. The encoded protein catalyzes this reaction at low concentrations of magnesium chloride while the other enzyme has high activity at high concentrations of magnesium chloride. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],</p>