

DSG4 Polyclonal Antibody

Cat No: HR1AP10698

For research use only

Overview

Product Name	DSG4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	114kDa
GeneID?Human?	DSG4 CDHF13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>desmoglein 4(DSG4) Homo sapiens This gene encodes a member of the desmoglein subgroup of desmosomal cadherins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein is a transmembrane component of desmosomes and may play a role in cell-cell adhesion in epithelial cells. Mutations in the gene are associated with localized autosomal recessive hypotrichosis and monilethrix, characterized by impaired hair growth. [provided by RefSeq, May 2016],</p>