

DHRS9 Polyclonal Antibody

Cat No: HR1AP10696

For research use only

Overview

Product Name	DHRS9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	DHRS9 UNQ835/PRO1773
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dehydrogenase/reductase 9(DHRS9) Homo sapiens This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. This protein demonstrates oxidoreductase activity toward hydroxysteroids and is able to convert 3-alpha-tetrahydroprogesterone to dihydroxyprogesterone and 3-alpha-androstenediol to dihydroxyprogesterone in the cytoplasm, and may additionally function as a transcriptional repressor in the nucleus. [provided by RefSeq, Jan 2014],</p>