

DHB14 Polyclonal Antibody

Cat No: HR1AP10694

For research use only

Overview

Product Name	DHB14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	HSD17B14 DHRS10 SDR3 UNQ502/PRO474
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hydroxysteroid 17-beta dehydrogenase 14(HSD17B14) Homo sapiens 17-beta-hydroxysteroid dehydrogenases, such as HSD17B14, are primarily involved in metabolism of steroids at the C17 position and also of other substrates, such as fatty acids, prostaglandins, and xenobiotics (Lukacik et al., 2007 [PubMed 17067289]).[supplied by OMIM, Jun 2009],