

DHRS7 Polyclonal Antibody

Cat No: HR1AP10695

For research use only

Overview

Product Name	DHRS7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	DHRS7 DHRS7A RETSDR4 CGI-86 UNQ285/PRO3448
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	dehydrogenase/reductase 7(DHRS7) Homo sapiens This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics. [provided by RefSeq, Apr 2016],