

DOC11 Polyclonal Antibody

Cat No: HR1AP10692

For research use only

Overview

Product Name	DOC11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	228kDa
GeneID?Human?	DOCK11 ZIZ2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:The DHR-2 domain is necessary for the GEF activity.,function:Guanine nucleotide-exchange factor (GEF) that activates CDC42 by exchanging bound GDP for free GTP.,miscellaneous:'Zizim' means 'spike' in Hebrew.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the DOCK family.,similarity:Contains 1 DHR-1 (CZH-1) domain.,similarity:Contains 1 DHR-2 (CZH-2) domain.,similarity:Contains 1 PH domain.,subunit:Interacts with CDC42.,</p>