

## 5NTC Polyclonal Antibody

Cat No: HR1AP10691

For research use only

### Overview

Product Name	5NTC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	NT5C2 NT5B NT5CP PNT5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	5'-nucleotidase, cytosolic II(NT5C2) Homo sapiens This gene encodes a hydrolase that serves as an important role in cellular purine metabolism by acting primarily on inosine 5'-monophosphate and other purine nucleotides. [provided by RefSeq, Oct 2011],