

## CNDP2 Polyclonal Antibody

Cat No: HR1AP10688

For research use only

### Overview

Product Name	CNDP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	CNDP2 CN2 CPGL PEPA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CNDP dipeptidase 2 (metallopeptidase M20 family)(CNDP2) Homo sapiens CNDP2, also known as tissue carnosinase and peptidase A (EC 3.4.13.18), is a nonspecific dipeptidase rather than a selective carnosinase (Teufel et al., 2003 [PubMed 12473676]).[supplied by OMIM, Mar 2008],