

## DC1I2 Polyclonal Antibody

Cat No: HR1AP10686

For research use only

### Overview

Product Name	DC1I2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	DYNC1I2 DNCI2 DNCIC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	dynein cytoplasmic 1 intermediate chain 2(DYNC1I2) Homo sapiens This gene encodes a member of the dynein intermediate chain family. The encoded protein is a non-catalytic component of the cytoplasmic dynein 1 complex, which acts as a retrograde microtubule motor to transport organelles and vesicles. A pseudogene of this gene is located on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012],