

CNGA3 Polyclonal Antibody

Cat No: HR1AP10680

For research use only

Overview

Product Name	CNGA3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	CNGA3 CNCG3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cyclic nucleotide gated channel alpha 3(CNGA3) Homo sapiens This gene encodes a member of the cyclic nucleotide-gated cation channel protein family which is required for normal vision and olfactory signal transduction. Mutations in this gene are associated with achromatopsia (rod monochromacy) and color blindness. Two alternatively spliced transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008],