

CSN7A Polyclonal Antibody

Cat No: HR1AP10678

For research use only

Overview

Product Name	CSN7A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	COPS7A CSN7A DERP10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	COP9 signalosome subunit 7A(COPS7A) Homo sapiens This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described. [provided by RefSeq, Mar 2014],