

CNTN3 Polyclonal Antibody

Cat No: HR1AP10675

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | CNTN3 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 113kDa |
| GeneID?Human? | CNTN3 KIAA1496 PANG |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting activity.,similarity:Belongs to the immunoglobulin superfamily. Contactin family.,similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 6 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:In brain, it is expressed in frontal lobe, occipital lobe, cerebellum and amygdala., |