

CNDD3 Polyclonal Antibody

Cat No: HR1AP10674

For research use only

Overview

Product Name	CNDD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	164kDa
GeneID?Human?	NCAPD3 CAPD3 KIAA0056
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>non-SMC condensin II complex subunit D3(NCAPD3) Homo sapiens Condensin complexes I and II play essential roles in mitotic chromosome assembly and segregation. Both condensins contain 2 invariant structural maintenance of chromosome (SMC) subunits, SMC2 (MIM 605576) and SMC4 (MIM 605575), but they contain different sets of non-SMC subunits. NCAPD3 is 1 of 3 non-SMC subunits that define condensin II (Ono et al., 2003 [PubMed 14532007]).[supplied by OMIM, Mar 2008].</p>