

C1T9B Polyclonal Antibody

Cat No: HR1AP10670

For research use only

Overview

Product Name	C1T9B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	C1QTNF9B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 C1q domain.,similarity:Contains 3 collagen-like domains.,