

CCDC8 Polyclonal Antibody

Cat No: HR1AP10668

For research use only

Overview

Product Name	CCDC8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	CCDC8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>coiled-coil domain containing 8(CCDC8) Homo sapiens This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a cofactor required for p53-mediated apoptosis following DNA damage, and may also play a role in growth through interactions with the cytoskeletal adaptor protein obscurin-like 1. Mutations in this gene are a cause of 3M syndrome-3 (3M3). [provided by RefSeq, Dec 2011],</p>