

MYCPP Polyclonal Antibody

Cat No: HR1AP10665

For research use only

Overview

Product Name	MYCPP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	204kDa
GeneID?Human?	DENND4A IRLB MYCPBP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>DENN domain containing 4A(DENND4A) Homo sapiens This gene encodes a DENN domain-containing protein that may function as a guanine nucleotide exchange factor that specifically activates ras-related protein Rab-10. This protein also contains a interferon stimulated response element-binding domain and may be involved in regulating the v-myc avian myelocytomatosis viral (MYC) oncogene. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 8. [provided by RefSeq, Mar 2016],</p>