

## CPSF6 Polyclonal Antibody

Cat No: HR1AP10663

For research use only

### Overview

Product Name	CPSF6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	CPSF6 CFIM68
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cleavage and polyadenylation specific factor 6(CPSF6) Homo sapiens The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. The cleavage factor complex is composed of four polypeptides. This gene encodes the 68kD subunit. It has a domain organization reminiscent of spliceosomal proteins. [provided by RefSeq, Jul 2008].</p>