



KIF4B Polyclonal Antibody

Cat No: HR1AP10653

For research use only

Overview

Product Name	KIF4B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	135kDa
GeneID?Human?	KIF4B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	kinesin family member 4B(KIF4B) Homo sapiens This gene is an intronless retrocopy of kinesin family member 4A. The protein encoded by this gene is a microtubule-based motor protein that plays vital roles in anaphase spindle dynamics and cytokinesis. [provided by RefSeq, Jul 2016].