

CHD5 Polyclonal Antibody

Cat No: HR1AP10649

For research use only

Overview

Product Name	CHD5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	214kDa
GeneID?Human?	CHD5 KIAA0444
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromodomain helicase DNA binding protein 5(CHD5) Homo sapiens This gene encodes a member of the chromodomain helicase DNA-binding protein family. Members of this family are characterized by a chromodomain, a helicase ATP-binding domain and an additional functional domain. This gene encodes a neuron-specific protein that may function in chromatin remodeling and gene transcription. This gene is a potential tumor suppressor gene that may play a role in the development of neuroblastoma. [provided by RefSeq, Feb 2012],</p>