

CHD2 Polyclonal Antibody

Cat No: HR1AP10648

For research use only

Overview

Product Name	CHD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	201kDa
GeneID?Human?	CHD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromodomain helicase DNA binding protein 2(CHD2) Homo sapiens The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>