

CHSS3 Polyclonal Antibody

Cat No: HR1AP10646

For research use only

Overview

Product Name	CHSS3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	97kDa
GeneID?Human?	CHSY3 CHSY2 CSS3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	chondroitin sulfate synthase 3(CHSY3) Homo sapiens CSS3 is a glycosyltransferase that has both glucuronyltransferase and N-acetylgalactosaminyltransferase activities (Yada et al., 2003 [PubMed 12907687]).[supplied by OMIM, Mar 2008],