

CLCN2 Polyclonal Antibody

Cat No: HR1AP10642

For research use only

Overview

Product Name	CLCN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	98kDa
GeneID?Human?	CLCN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	chloride voltage-gated channel 2(CLCN2) Homo sapiens This gene encodes a voltage-gated chloride channel. The encoded protein is a transmembrane protein that maintains chloride ion homeostasis in various cells. Defects in this gene may be a cause of certain epilepsies. Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2012],