

## PDE3B Polyclonal Antibody

Cat No: HR1AP10640

For research use only

### Overview

Product Name	PDE3B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	PDE3B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:Nucleoside 3',5'-cyclic phosphate + H(2)O = nucleoside 5'-phosphate.,enzyme regulation:Inhibited by cGMP.,function:May play a role in fat metabolism.,similarity:Belongs to the cyclic nucleotide phosphodiesterase family.,tissue specificity:Abundant in adipose tissues.,