

CDR1 Polyclonal Antibody

Cat No: HR1AP10637

For research use only

Overview

Product Name	CDR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	CDR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cerebellar degeneration related protein 1(CDR1) Homo sapiens Autoantibodies directed against the protein encoded by this intronless gene have been found in some patients with paraneoplastic cerebellar degeneration. The encoded protein contains several hexapeptide repeats. [provided by RefSeq, Jan 2010],