

CEP68 Polyclonal Antibody

Cat No: HR1AP10633

For research use only

Overview

Product Name	CEP68 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	CEP68 KIAA0582
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CEP68 (Centrosomal Protein 68) is a Protein Coding gene. Diseases associated with CEP68 include Retinitis Pigmentosa 28. Involved in maintenance of centrosome cohesion, probably as part of a linker structure which prevents centrosome splitting (PubMed: 18042621). Required for localization of CDK5RAP2 to the centrosome during interphase.