

CEP72 Polyclonal Antibody

Cat No: HR1AP10634

For research use only

Overview

Product Name	CEP72 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	71kDa
GeneID?Human?	CEP72 KIAA1519
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centrosomal protein 72(CEP72) Homo sapiens The product of this gene is a member of the leucine-rich-repeat (LRR) superfamily of proteins. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. [provided by RefSeq, Jul 2008],