

HAUS2 Polyclonal Antibody

Cat No: HR1AP10631

For research use only

Overview

Product Name	HAUS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	HAUS2 C15orf25 CEP27
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	HAUS augmin like complex subunit 2(HAUS2) Homo sapiens The protein encoded by this gene is a subunit of the augmin complex. The augmin complex plays a role in microtubule attachment to the kinetochore and central spindle formation. [provided by RefSeq, Apr 2016],