

CEP63 Polyclonal Antibody

Cat No: HR1AP10632

For research use only

Overview

Product Name	CEP63 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	CEP63
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centrosomal protein 63(CEP63) Homo sapiens This gene encodes a protein with six coiled-coil domains. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. Several alternatively spliced transcript variants have been found, but their biological validity has not been determined. [provided by RefSeq, Jul 2008],