

ACAP2 Polyclonal Antibody

Cat No: HR1AP10630

For research use only

Overview

Product Name	ACAP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	ACAP2 CENTB2 KIAA0041
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	enzyme regulation:GAP activity stimulated by phosphatidylinositol 4,5-bisphosphate (PIP2) and phosphatidic acid.,function:GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 Arf-GAP domain.,similarity:Contains 1 BAR domain.,similarity:Contains 1 PH domain.,similarity:Contains 3 ANK repeats.,tissue specificity:Widely expressed. Highest level in lung.,