

BORG5 Polyclonal Antibody

Cat No: HR1AP10628

For research use only

Overview

Product Name	BORG5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	CDC42EP1 BORG5 CEP1 MSE55
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CDC42 effector protein 1(CDC42EP1) Homo sapiens CDC42 is a member of the Rho GTPase family that regulates multiple cellular activities, including actin polymerization. The protein encoded by this gene is a CDC42 binding protein that mediates actin cytoskeleton reorganization at the plasma membrane. This protein is secreted and is primarily found in bone marrow. [provided by RefSeq, Jul 2008],