

CTR3 Polyclonal Antibody

Cat No: HR1AP10627

For research use only

Overview

Product Name	CTR3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	SLC7A3 ATRC3 CAT3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 7 member 3(SLC7A3) Homo sapiens This gene encodes a member of the solute carrier family 7. The encoded protein is a sodium-independent cationic amino acid transporter. Alternate splicing results in multiple transcripts that encoded the same protein.[provided by RefSeq, May 2010],