

CTSRB Polyclonal Antibody

Cat No: HR1AP10624

For research use only

Overview

Product Name	CTSRB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	CATSPERB C14orf161
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probably involved in sperm cell hyperactivation via its association with CATSPER1. Sperm cell hyperactivation is needed for sperm motility which is essential late in the preparation of sperm for fertilization.,subunit:Component of a complex, at least composed of CATSPER1, HSPA1B and CATSPERB.,