

CTSR1 Polyclonal Antibody

Cat No: HR1AP10621

For research use only

Overview

Product Name	CTSR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	CATSPER1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cation channel sperm associated 1(CATSPER1) Homo sapiens Calcium ions play a primary role in the regulation of sperm motility. This gene belongs to a family of putative cation channels that are specific to spermatozoa and localize to the flagellum. The protein family features a single repeat with six membrane-spanning segments and a predicted calcium-selective pore region. [provided by RefSeq, Jul 2008],</p>