

CEAM7 Polyclonal Antibody

Cat No: HR1AP10618

For research use only

Overview

Product Name	CEAM7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	CEACAM7 CGM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carcinoembryonic antigen related cell adhesion molecule 7(CEACAM7) Homo sapiens This gene encodes a cell surface glycoprotein and member of the carcinoembryonic antigen (CEA) family of proteins. Expression of this gene may be downregulated in colon and rectal cancer. Additionally, lower expression levels of this gene may be predictive of rectal cancer recurrence. This gene is present in a CEA family gene cluster on chromosome 19. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2015],</p>