

CEA20 Polyclonal Antibody

Cat No: HR1AP10617

For research use only

Overview

Product Name	CEA20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	CEACAM20 UNQ9366/PRO34155
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the immunoglobulin superfamily. CEA family.,similarity:Contains 4 Ig-like C2-type (immunoglobulin-like) domains.,