

CEA19 Polyclonal Antibody

Cat No: HR1AP10616

For research use only

Overview

Product Name	CEA19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	CEACAM19 CEAL1 UNQ2973/PRO7436
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	induction:Down-regulated by dexamethasone (in vitro).,similarity:Belongs to the immunoglobulin superfamily. CEA family.,tissue specificity:Ubiquitous with highest expression in prostate, uterus, fetal brain, mammary gland, adrenal gland, skeletal muscle, small intestine, and kidney, and lower expression in lung, cerebellum, testis, liver, pancreas, bone marrow and ovary.,