



## CBPO Polyclonal Antibody

Cat No: HR1AP10614

For research use only

### Overview

Product Name	CBPO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	CPO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	carboxypeptidase O(CPO) Homo sapiens This gene is a member of the metallo-carboxypeptidase gene family. [provided by RefSeq, Jan 2011],