

CBPD Polyclonal Antibody

Cat No: HR1AP10613

For research use only

Overview

Product Name	CBPD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	151kDa
GeneID?Human?	CPD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carboxypeptidase D(CPD) Homo sapiens The metallo-carboxypeptidase family of enzymes is divided into 2 subfamilies based on sequence similarities. The pancreatic carboxypeptidase-like and the regulatory B-type carboxypeptidase subfamilies. Carboxypeptidase D has been identified as a regulatory B-type carboxypeptidase. CPD is a homolog of duck gp180, a hepatitis B virus-binding protein. Transcript variants utilizing alternative polyadenylation signals exist for this gene. [provided by RefSeq, Jul 2008],</p>