

# CBPB1 Polyclonal Antibody

Cat No: HR1AP10612

For research use only

## Overview

Product Name	CBPB1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	CPB1 CPB PCPB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carboxypeptidase B1(CPB1) Homo sapiens Three different procarboxypeptidases A and two different procarboxypeptidases B have been isolated. The B1 and B2 forms differ from each other mainly in isoelectric point. Carboxypeptidase B1 is a highly tissue-specific protein and is a useful serum marker for acute pancreatitis and dysfunction of pancreatic transplants. It is not elevated in pancreatic carcinoma. [provided by RefSeq, Jul 2008],</p>