

CBPA6 Polyclonal Antibody

Cat No: HR1AP10611

For research use only

Overview

Product Name	CBPA6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	CPA6 CPAH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carboxypeptidase A6(CPA6) Homo sapiens The gene encodes a member of the peptidase M14 family of metalloproteinases. The encoded preproprotein is proteolytically processed to generate the mature enzyme, which catalyzes the release of large hydrophobic C-terminal amino acids. This enzyme has functions ranging from digestion of food to selective biosynthesis of neuroendocrine peptides. Mutations in this gene may be linked to epilepsy and febrile seizures, and a translocation t(6;8)(q26;q13) involving this gene has been associated with Duane retraction syndrome. [provided by RefSeq, May 2016].</p>